

Date 8/19/2013

Objective Rd 2 PCR amplification of triplet repeat region of FMR1 using the Rd 1 PCR pdt as template:

- 1:10 and 1:50 dilutions of Rd 1 PCR pdt advised
- Cycling conditions of Rd 1 modified as:
 - Three step amplification process
 - Denat temp set to 92degC/ 30secs
 - Annealing temp grad* 50-70deg/ 1:30mins
 - Extension temp. set to 72degC/ 2:30mins

***: PCR conditions as follows - 90/70/72 has worked before
(Oct-Nov 2012) with taq for rd 2 PCR**

In fact ann'g at temps 70, 68, 66 and 62 have all worked.

Description *Template used:* Rd 1 PCR pdt to a final conc of 0.1X or 0.02X.
Primers: FT for and FT rev
Polymerase: Native Taq Polymerase (Life Technologies)
Solvents: Rxns in 0.5M NMP + 2.2M Betaine

Reaction Composition:

Component	0.5M NMP + 2.2M Betaine rxns Final Conc/ amt
Taq Buffer	1 X
MgCl ₂	1.5 mM
dNTPs	0.4 mM
FT-For/ SFS-for	0.4 uM
FT-Rev/ SFS-rev	0.4 uM
Rd 1 PCR rxn	0.1X or 0.02X
Taq Polymerase	2.5 U
Water	
NMP	0.5 M
Betaine	2.2 M
total	25ul

Cycling Conditions:


- 1 92degC/ 30secs
- 2 50-70degC/ 1:30mins
- 3 72degC/ 2:30mins
- 4 GOTO 2 40 times

Gel Electrophoresis:

10% precast polyacrylamide gels (Life Technologies)
 8ul of PCR rxn + 2ul of (5X) sample buffer loaded
 Molecular Ladders:
 100bp ladder (NEB) 5ul
 50bp ladder (NEB) 5ul
 20bp ladder (BioRad) 3ul

Gel Pictures

FT primers - NMP + Betaine-no hot start(20sec exposure)



Gel Pictures

FT primers - NMP + Betaine-with hot start (20sec exposure)



Comments None of the reactions have worked.

When attempting to reproduce Ram's work in Oct-Nov 2012, Rd 2 PCR using 0.1X, 0.02X, 0.01X of Rd 1 PCR with Taq with 90-70-72 has worked to produce a PCR ptd.

See power point file: FMR-FT PCR Amplification Oct-Nov 2012 121012 (slide 11)

Hence these conditions were tested again. However, the PCR did not work.

Round 2 PCR: NA20230 0.1-0.002X Round 1 reaction mix with native Taq Polymerase in 0.5M NMP +2.2M Betaine

