

Date 8/22/2013

Objective Rd 2 PCR amplification of triplet repeat region of FMR1 using the Rd 1 PCR pdt as template:

- 1:10 and 1:50 dilutions of Rd 1 PCR pdt advised
- Cycling conditions of Rd 1 modified as:
 - Three step amplification process
 - Denat temp set to 92degC/ 30secs
 - Annealing temp set to 70deg/ 1:30mins
 - Extension temp. set to 72degC/ 2:30mins

Description *Template used:* Rd 1 PCR pdt 0.1X, 0.02X, 0.01X
Primers: FT for and FT rev
Polymerase: Deep Vent (NEB)
Solvents: No Solvt, Rxns in 0.5M NMP, Rxns in 2.2M Betaine

Reaction Composition:

Component	No Solvt Final Conc/	0.5M NMP Final Conc/	2.2M Betaine Final Conc/
ThermoPol Buffer	1 X	1 X	1 X
dNTPs	0.4 mM	0.4 mM	0.4 mM
FT-For	0.4 uM	0.4 uM	0.4 uM
FT-Rev	0.4 uM	0.4 uM	0.4 uM
Rd 1 PCR pdt*	0.1 X	0.1 X	0.1 X
DV Polymerase	0.25 U	0.25 U	0.25 U
Water			
NMP		0.5 M	
Betaine			2.2 M
total	25.00	25.00	25.00

Cycling Conditions:


- 1 98deg/ 4mins
PAUSE-Add DV-RESUME
- 2 92degC/ 30secs
- 3 70degC/ 1:30mins
- 4 72degC/ 2:30mins
- 5 GOTO 2 40 times

Gel Electrophoresis:

10% precast polyacrylamide gels (Life Technologies)
 8ul of PCR rxn + 2ul of (5X) sample buffer loaded
 Molecular Ladders:
 100bp ladder (NEB) 5ul
 50bp ladder (NEB) 5ul
 20bp ladder (BioRad) 3ul

Gel Pictures

FT primers - No Solvt-NMP-Betaine-w/ hot start (20sec exposure)



Comments

Rd 2 PCR has worked better previously with 0.01X Rd 1 Pdt (using DV) in NMP or betaine
See power point file: FMR-FT PCR Amplification Oct-Nov 2012 121012 (slide 18)
It has not worked as well in the current expt, however, the only band seen in the target band.
This needs to be repeated to confirm.

Round 2 PCR (NA20230): Use of Deep Vent

- denaturation temperature gradient 90°-98°C with annealing set at 70°C and extension set at 72°C
- testing different solvent combinations: 0.5M NMP, 2.2M Betaine, 0.5M NMP + 2.2M Betaine,

