

Date 2/8/2013

Objective PCR amplification of triplet repeat region of FMR1 with a modified program (as below):
 - Without hot start
 - Three step amplification process
 - Denat temp grad 75-95deg
 - Extension and annealing temp.s and times (that have worked previously
 ie 54deg/ 1:30mins annealing-72deg/ 2:30mins extension)

Description *Template used:* 100ng genomic DNA NA20230, (53/54 rpt sample from Coriell)
 This is the newly ordered (Jan 2013) DNA: all old DNA is over
Primers: FT for and FT rev
Polymerase: Deep Vent (New England Biolabs)
Solvents: No Solvt rxns
 Rxns in 1M TMSO

<i>Reaction Composition:</i>	No solvt rxns		1M TMSO rxns	
	Component	Final Conc/amt	Final Conc/amt	
	Thermopol Buffer	1 X	1 X	
	dNTPs	0.4 mM	0.4 mM	
	FT-For/ SFS-for	0.4 uM	0.4 uM	
	FT-Rev/ SFS-rev	0.4 uM	0.4 uM	
	gDNA NA203230	100 ng	100 ng	
	Deep Vent	0.25 U	0.25 U	
	Water			
	TMSO		1 M	
	total	25ul	25ul	

Cycling Conditions:


- 1 75-95degC/ 30secs
- 2 54degC/ 1:30mins
- 3 72degC/ 2:30mins
- 4 GOTO 2 40 times

Gradient steps: 95, 94, 91.7, 87.7, 82.9, 79.1, 76.4, 75degC


Gel Electrophoresis: 10% precast polyacrylamide gels (Life Technologies)
 8ul of PCR rxn + 2ul of (5X) sample buffer loaded
 Molecular Ladders: 5ul loaded
 100bp ladder (NEB)
 50bp ladder (NEB)

Gel Pictures

Lanes are marked according to the Denat temp of the particular rxn
 Expected product size is 383bp
FT-DV- No Solvt 20sec exposure



FT-DV-1M TMSO 20sec exposure



Comments

No bands of expected size.
 Based on Taq results, expt will be repeated with new primer stks as well.