

Date 1/29/2013

Objective PCR amplification of triplet repeat region of FMR1 with a modified program (as below):
 - Without hot start
 - Denaturation temp gradient 75-95degC
 - Three step amplification process
 - Extension and annealing temp.s and times (that have worked previously
 ie 54deg/ 1:30mins annealing-72deg/ 2:30mins extension)

Description *Template used:* 100ng genomic DNA NA20230, (53/54 rpt sample from Coriell)
Primers: SFS for and SFS rev
 FT for and FT rev
Polymerase: Native Taq Polymerase (Life Technologies)
Solvents: No Solvt rxns
 Rxns in 0.5M NMP + 2.2M Betaine

Reaction Composition:

Component	No Solvt Rxns Final Conc/ amt	0.5M NMP + 2.2M Betaine rxns Final Conc/ amt
Taq Buffer	1 X	1 X
MgCl2	1.5 mM	1.5 mM
dNTPs	0.4 mM	0.4 mM
FT-For/ SFS-for	0.4 uM	0.4 uM
FT-Rev/ SFS-rev	0.4 uM	0.4 uM
gDNA NA203230	100 ng	100 ng
Taq Polymerase	2.5 U	2.5 U
Water		
NMP		0.5 M
Betaine		2.2 M
total	25ul	25ul

Cycling Conditions:

- 1 75-95degC/ 30secs
- 2 54degC/ 1:30mins
- 3 72degC/ 2:30mins
- 4 GOTO 2 40 times

Gradient steps: 95, 94, 91.7, 87.7, 82.9, 79.1, 76.4, 75degC

Gel Electrophoresis: 10% precast polyacrylamide gels (Life Technologies)
 8ul of PCR rxn + 2ul of (5X) sample buffer loaded
 Molecular Ladders: 5ul loaded
 100bp ladder (NEB)
 50bp ladder (NEB)

Gel Pictures

All pictures are 20sec exposures

Expected product size is 383 bp

Lanes are marked according to the denaturation temp of the particular rxn.


See **012913 NA20230 rePCRd FT-SFS PMC Program 2 denat 75-95**

for the gel pics of SFS PCRs


This expt originally done on 1/29/13 has been put here again (only FT primer data)

for purposes of completeness

FT primers - No solvt



FT primers - NMP + Betaine



- Comments** **Only rxns with FT primers in 0.5M NMP + 2.2M Betaine show expected bands:**
- Reactions without hot start work at a denat temp of 95, 94, 91.7degC giving the expected pdt of 383bp.
 - Lower denat temps do not work.
 - However, even in the positive rxns, there are multiple other bands (of both smaller and higher sizes).
 - Hotstart 1 and Hotstart 2 were set up from the same mastermix stored O/N (2.5U taq in N+B solvt O/N) at 4degC. However, in these rxns more taq was added after the hot start and before cycling
 - These also show the 383bp band but show many more extra bands.
 - The hot start rxns were PCR'd as follows:
 - denat 90deg/ 30secs
 - annealing 54deg/ 1:30mins
 - extension 72deg/ 2:30mins
 - 40 cycles
- All other test conditions do not produce any bands of expected size. Even for these corresponding rxns were set up with hot start (as described above). However, even the hot start rxns do not produce any bands of expected size. The 79.1deg denat rxn with FT primers (no solvt) shows a band in the expected region but is more likely smaller than 383bp. Not sure.
- Conclusion** 92degC/ 30secs is going to be the preferred denat temp (without hot start)