


SIRT3 Purification

1. SIRT3 plasmids were transformed to BL21 cell.
2. Plate the above chemical transformation product on LB argar with Ampicillin (1:1000) at 37oC overnight.
3. Single colonies were detected and were cultured firstly in 500 ul LB broth/Amp till OD to 0.6.
4. Transfer 20 ul of above culture to 5mL LB broth/Amp overnight at 37oC with 250 rpm
5. Scale up for 100ml of cell culture.
6. Add IPTG for final concentration (2mM) at 37oC 250 rpm for 3 hours.
7. Centrifuge to get the cell pellet and stored at -80oC till use.
8. Add lysis buffer + bugbuster and sonication.
9. Centrifuge to obtain supernatant

SIRT3 (102-399)
UniProt ID: Q9NTG7
Length: 298
Mass (Da): 32,700




Under current condition, protein expression level is good.

10. AKTA first Ni-column purification


- Equilibration 4.5 ml
- Wash 20 column volume (20 ml, 5ml /fraction)
- Elution 10 cv (10 ml, 0.5 ml/fraction)
- Waste

11. Run protein gel to check

1. SIRT3 protein expression level
2. Fraction content
 - First 2 wash (10 ml) are enough
 - First 7 elution (3.5ml) are good
3. Modify AKTA program
 - Reduce wash volume to 10 CV
 - Reduce elution to 5 CV



- Raw sequencing data for SIRT1 and SIRT2 have been obtained from OriGene.
- The circular maps have been drawn.
- The restriction sites have been obtained
 - SIRT1: Kpn1 and Not1
 - SIRT2: Ecor1 and Not1
- Double digestion was performed for SIRT1 and SIRT2



- Since different restriction enzymes have different working buffer, the double digest Finder from NEB web site was used to find right buffer.



Double Digest Finder

SIRT1

Use this tool to guide your reaction buffer selection when setting up double digests, a common timesaving procedure. Choosing the right buffers will help you to avoid star activity and loss of product. [Learn more information about double digests](#), including how to set up a reaction.

Select 1st enzyme

 Select 2nd enzyme

Enzyme	Cat #	Temp	Supplied NEBuffer	Supplements	% Activity in NEBuffer			
					SAM	1.1	2.1	3.1
KpnI 	R0142	37°C	NEBuffer 1.1	no	100	75	10	50*
NotI 	R0189	37°C	NEBuffer 3.1	no	10	50	100	25



Double Digest Finder

SIRT2

Use this tool to guide your reaction buffer selection when setting up double digests, a common timesaving procedure. Choosing the right buffers will help you to avoid star activity and loss of product. [Learn more information about double digests](#), including how to set up a reaction.


Select 1st enzyme

 Select 2nd enzyme

Enzyme	Cat #	Temp	Supplied NEBuffer	Supplements	% Activity in NEBuffer			
					SAM	1.1	2.1	3.1
EcoRI 	R0101	37°C	NEBuffer EcoRI	no	25	100*	50	50*
NotI 	R0189	37°C	NEBuffer 3.1	no	10	50	100	25

- SIRT1 NEBuffer 2.1
- SIRT2 NEBuffer 3.1
- New protein ladder was ordered.
- NEBuffer 2.1 was ordered.

SIRT1



SIRT2

